

# Large-Volume | Smart Device



[www.simatec-medical.com](http://www.simatec-medical.com)


 **SWISS MADE**



# Large-Volume | Smart Device



## Device specification

Drive	Gas producing dry cell
Volume	20ml / 50ml / 100ml*
Viscosity	1..>100 cP
Flow rate	0.5-30ml/h (fix or adjustable)* Starts/Stops at any time on demand
Accuracy	< +/-5%
Weight	30g, drive system
Handling	Pre-filled or patient-filled
Temperature	-20°C..+55°C
Storage time	>3y of the gas producing cell
Primary container	Polymer / glass
Drive benefits	Silent; pulsation free
Pressure	25bar
Waste disposal	Battery recycling 
*	Other specifications on request

[www.simatec-medical.com](http://www.simatec-medical.com)